

## Determination of Selected Rice Herbicides in Groundwater

### 1. Scope:

This section method (SM) is applicable to the analysis of Bensulfuron-methyl, Bispyribac sodium, Clomazone, Halosulfuron methyl, MCPA, Molinate, Orthosulfamuron, Penoxsulam, Propanil, Propiconazole, Thiobencarb and Triclopyr in groundwater. It is followed by all authorized EMON personnel.

### 2. Principle:

Groundwater sample (500 mL) is acidified with 3 N Hydrochloric acid. A solid phase extraction HLB cartridge is used to retain the selected rice herbicides from acidified water samples. The analytes are eluted with acetonitrile/methanol (50:50) solution. The eluant is concentrated and analyzed by Liquid Chromatography coupled to a Linear Ion Trap Quadrupole LC/MS/MS. The reporting limit is 0.05 ppb for all compounds.

### 3. Safety:

All general laboratory safety rules for sample preparation and analysis shall be followed.

### 4. Interferences:

There were no matrix interferences that caused quantitative problems during method development and validation.

### 5. Apparatus and Equipment:

- 5.1 Nitrogen Evaporator (Meyer N-EVAP Organomation Model #112 or equivalent)
- 5.2 Balance (Mettler PC 4400 or equivalent)
- 5.3 Vortex-vibrating mixer
- 5.4 Solid phase extraction manifold, Supelco Visiprep TM24 or equivalent
- 5.5 Solid phase extraction manifold accessories: vacuum source, vacuum chamber, vacuum controller, cartridge fittings (tube adapters) and connectors, sample delivery tubing with stainless steel weight, sample collection tubes and rack.
- 5.6 pH meter Mettler-Toledo or equivalent
- 5.7 Ultra High-Performance Liquid Chromatography (UHPLC) equipped with a linear ion trap quadrupole (MS/MS)

6. Reagents and Supplies:

- |      |   |                  |
|------|---|------------------|
| 6.1  | Bensulfuron-methyl  | CAS# 83055-99-6  |
| 6.2  | Bispyribac sodium   | CAS# 125401-92-5 |
| 6.3  | Clomazone   | CAS# 81777-89-1  |
| 6.4  | Halosulfuron-methyl   | CAS#100784-20-1  |
| 6.5  | MCPA  | CAS#94-74-6      |
| 6.6  | Molinate  | CAS#2212-67-1    |
| 6.7  | Orthosulfamuron   | CAS#213464-77-8  |
| 6.8  | Penoxsulam  | CAS#219714-96-2  |
| 6.9  | Propanil  | CAS#709-98-8     |
| 6.10 | Propiconazole   | CAS#60207-90-1   |
| 6.11 | Thiobencarb   | CAS#28249-77-6   |
| 6.12 | Triclopyr   | CAS#55335-06-3   |
| 6.13 | 2,4,5 T   | CAS#93-76-5      |
| 6.14 | Simazine-d10  | CAS#220621-39-6  |
| 6.15 | Diclofenac-d4 (IS)  | CAS#153466-65-0  |
| 6.16 | Methanol, MS grade, Fisher Optima LC/MS or equivalent   |                  |
| 6.17 | Water, MS grade, Fisher Optima LC/MS or equivalent  |                  |
| 6.18 | Acetonitrile, Fisher Optima or equivalent   |                  |
| 6.19 | Hydrochloric acid 3 N   |                  |
| 6.20 | Pipettes; air cushioned and positive displacement, various volumes and types, Eppendorf or equivalent |                  |
| 6.21 | Disposable Pasteur pipettes and other laboratory ware as needed                                       |                  |
| 6.22 | Solid phase extraction cartridges: Waters Oasis® HLB 6 cc (200 mg),                                   |                  |
| 6.23 | Graduated test tube, 15 mL (calibrated at 0.5mL with methanol)  |                  |
| 6.24 | LCMS Columns:<br>Analytical column: ACE Excel C18 2.0 µm, 2.1 x 100mm column or equivalent            |                  |
| 6.25 | Formic acid, HPLC grade   |                  |
| 6.26 | Aqueous Solution: Water with 0.1% Formic acid   |                  |
| 6.27 | Organic Solution: Methanol/Acetonitrile (50/50) with 0.1% Formic acid                                 |                  |

7. Standards Preparation:

- 7.1 Individual stock standards of 1.0mg/mL were obtained from the CDFA/CAC Standards Repository. Standards are ordered from a 17034 accredited supplier when possible.

The standards were diluted to 10 µg/mL with acetonitrile. A combination standard of 10 µg/mL was prepared from the individual mg/mL standards in acetonitrile. The combination standard was also used to dilute to the following concentrations: 0.05, 0.1, 0.2, 0.5, 1.0 and 2.0 µg/mL in acetonitrile which were later diluted 1:10 with clean background matrix for instrument calibration.

7.2 Keep all standards in the designated freezer for storage.

7.3 The expiration date of each standard is six months from the preparation date.

## 8. Sample Preservation and Storage:

Store all samples waiting for extraction in a designated refrigerator ( $4 \pm 3$  °C). Sample storage location and movement shall be recorded.

## 9. Test Sample Preparation:

### 9.1 Background Preparation

The Department of Pesticide Regulation (DPR) provided the groundwater for background to be used in method validation and QC.

### 9.2 Preparation of blank and spike

Matrix blank: Weigh out 500 g of background water and follow the test sample extraction procedure.

Matrix spike: Weigh out 500 g of background water. Spike a client requested amount of pesticide (0.1 – 0.2 ppb) into the background water, mix well and let it stand for one minute. Follow the test sample extraction procedure.

### 9.3 Test Sample Extraction

9.3.1 Remove samples from refrigerator and allow them to come to ambient temperature.

9.3.2 Weigh  $500 \pm 0.5$  g of water sample into a 600 mL beaker.

9.3.3 Add 0.1 µg 2,4,5 T / Simazine-d10 (100 µL of 1 ng/µL spiking solution) as a surrogate to each sample except blank.

- 9.3.4 Adjust pH to 2.0 – 2.5 with 3 N HCL.
- 9.3.5 Connect an HLB cartridge (200mg) to the vacuum manifold.
- 9.3.6 Condition the cartridge with ~5 mL of methanol at a flow rate ~ 8 mL/minutes followed by ~ 5 mL of acidified D.I. water (~ pH 2) by applying vacuum.
- 9.3.7 Turn off the vacuum when the D.I. water has just passed through the cartridges. Refill HLB cartridges with acidified D.I. water. Attach the sample delivery tubes to the cartridge and place weighted tube ends into water sample.
- 9.3.8 Allow the sample to pass through the conditioned cartridges by applying vacuum. Adjust the flow rate to ~ 8 mL/minute
- 9.3.9 After all the water sample has passed through the cartridges, increase the vacuum to ~ 20 psi for 10 minutes. Detach the sample delivery tube from HLB cartridge. Shake out any excess water in the cartridge reservoir and dry inside of cartridge with Kimwipe.
- 9.3.10 Place the graduated test tubes into the vacuum manifold.
- 9.3.11 Elute and collect all analytes with 5 mL of acetonitrile/methanol (50/50) solution and at a flow rate of ~8 mL/minutes. Concentrate the eluant to ~1.8 mL in a water bath at  $38 \pm 2$  °C under a gentle stream of nitrogen. Record the temperature into extraction sheet. Bring to a final volume of 2.0 mL with acetonitrile/methanol (50/50) solution.
- 9.3.13 Transfer the extract into an autosampler vial to be analyzed by ESI/LC/MS/MS.

## 10. Instrument Calibration:

- 10.1 A minimum of 3 levels of standards is required for linear curve.
- 10.2 The quadratic calibration standard curve consists of a minimum of five levels.

The recommended concentration levels of standards are 0.005, 0.01, 0.02, 0.05, 0.1, 0.2 µg/mL.

- 10.3 The calibration may be achieved using linear or quadratic regression with a correlation coefficient ( $r$ )  $\geq 0.995$  or ( $r^2$ )  $\geq 0.990$ .

## 11. Analysis:

### 11.1 Injection Scheme

The instrument may need to be conditioned with a few standards before running the following sequence of Standard Curve, Solvent, Matrix Blank, Matrix Spike, Test Samples and Standard Curve.

### 11.2 Linear Ion Trap Quadrupole LC/MS/MS Mass Spectrometer

#### 11.1.1 LC Instrument: Shimadzu LC30

Column: ACE Excel C18 2.0 µm, 2.1 x 100 mm column

Column Temperature: 40 °C

Mobile Phase: Gradient

Solvent 1: Aqueous Solution

Solvent 2: Organic Solution

| <u>Time (min)</u> | <u>Flow rate (mL/min)</u> | <u>Solvent 1</u> | <u>Solvent 2</u> |
|-------------------|---------------------------|------------------|------------------|
| 0.01              | 0.3                       | 65               | 35               |
| 2.0               | 0.3                       | 65               | 35               |
| 12.0              | 0.3                       | 15               | 85               |
| 15.0              | 0.3                       | 15               | 85               |
| 16.0              | 0.3                       | 65               | 35               |
| 17.0              | 0.3                       | 65               | 35               |

Injection Volume: 3.0 µL

#### 11.2.2 Mass Spectrometer and Operating Parameters

Model: ABSciex QTRAP 6500

Experiment 1

Ion ProbeType: Electrospray Ionization (ESI)

Ion Mode: Positive

Curtain Gas: 30  
Ion Spray Voltage: 4500  
Temp: 350  
Ion Source Gas 1 40  
Ion Source Gas 2 40  
Collision Gas: High

| Compound                        | RT    | Precursor Ion | Product Ion  | Declustering Potential | Collision Energy | Entrance Potential | Exit Potential |
|---------------------------------|-------|---------------|--------------|------------------------|------------------|--------------------|----------------|
| Bensulfuron-methyl              | 8.14  | 410.9         | <b>149</b>   | 66                     | 25               | 10                 | 18             |
|                                 |       | 410.9         | 182          | 66                     | 25               | 10                 | 24             |
| Bispyribac sodium               | 8.86  | 430.9         | <b>275</b>   | 51                     | 17               | 10                 | 14             |
|                                 |       | 430.9         | 413          | 51                     | 23               | 10                 | 20             |
| Clomazone                       | 7.95  | 240           | <b>125</b>   | 26                     | 23               | 10                 | 14             |
|                                 |       | 240           | 89.1         | 26                     | 65               | 10                 | 10             |
| Halosulfuron methyl             | 9.52  | 434.9         | <b>182</b>   | 61                     | 25               | 10                 | 22             |
|                                 |       | 434.9         | 138.9        | 61                     | 63               | 10                 | 22             |
| Molinate                        | 9.14  | 187.9         | <b>126</b>   | 41                     | 17               | 10                 | 16             |
|                                 |       | 187.9         | 55           | 41                     | 35               | 10                 | 14             |
| Orthosulfamuron                 | 7.38  | 424.9         | <b>199</b>   | 51                     | 17               | 10                 | 26             |
|                                 |       | 424.9         | 226.9        | 51                     | 19               | 10                 | 30             |
| Penoxsulam                      | 7.28  | 483.9         | <b>195</b>   | 80                     | 37               | 10                 | 26             |
|                                 |       | 483.9         | 443.9        | 80                     | 33               | 10                 | 20             |
| Propanil                        | 8.22  | 217.8         | <b>161.9</b> | 56                     | 21               | 10                 | 18             |
|                                 |       | 217.8         | 126.9        | 56                     | 33               | 10                 | 16             |
| Propiconazole                   | 11.06 | 341.9         | 158.9        | 56                     | 31               | 10                 | 20             |
|                                 |       | 341.9         | 69           | 56                     | 23               | 10                 | 12             |
| Thiobencarb                     | 11.51 | 257.9         | <b>124.9</b> | 41                     | 23               | 10                 | 16             |
|                                 |       | 257.9         | 89           | 41                     | 65               | 10                 | 14             |
| simazine-d10<br>(surrogate)     | 4.52  | 212.1         | <b>136.8</b> | 86                     | 27               | 10                 | 16             |
|                                 |       | 212.1         | 104.9        | 86                     | 35               | 10                 | 16             |
| Diclofenac-d4<br>(internal std) | 10.25 | 299.8         | <b>218.9</b> | 31                     | 27               | 10                 | 28             |
|                                 |       | 299.8         | 253.9        | 31                     | 19               | 10                 | 32             |

Experiment 2  
Ion ProbeType: Electropray Ionization (ESI)  
Ion Mode: Negative  
Curtain Gas: 30  
Ion Spray Voltage: -4500  
Temp: 350  
Ion Source Gas 1: 40  
Ion Source Gas 2: 40  
Collision Gas: High

| Compound                        | RT    | Precursor Ion | Product Ion  | Declustering Potential | Collision Energy | Entrance Potential | Exit Potential |
|---------------------------------|-------|---------------|--------------|------------------------|------------------|--------------------|----------------|
| MCPA                            | 7.54  | 199           | <b>141</b>   | -25                    | -18              | -10                | -9             |
|                                 |       | 199           | 155          | -25                    | -12              | -10                | -9             |
| Triclopyr                       | 8.13  | 254           | <b>195.9</b> | -10                    | -14              | -10                | -11            |
|                                 |       | 254           | 217.8        | -10                    | -8               | -10                | -13            |
| 2,4,5 T(surrogate)              | 8.93  | 253           | <b>194.8</b> | -10                    | -16              | -10                | -11            |
|                                 |       | 253           | 158.8        | -10                    | -38              | -10                | -11            |
| Diclofenac-d4<br>(internal std) | 10.25 | 298           | <b>254</b>   | -20                    | -16              | -10                | -15            |
|                                 |       | 298           | 217.1        | -20                    | -28              | -10                | -13            |

## 12. Quality Control:

### 12.1 Method Detection Limits (MDL)

Method Detection Limit (MDL) refers to the lowest concentration of the analyte that a method can detect reliably. To determine the MDL, 7 groundwater samples were spiked at 0.1ppb for the selected rice herbicides and processed through the entire method along with a blank. The standard deviation derived from the spiked sample recoveries was used to calculate the MDL using the following equation:

$$MDL = tS$$

Where t is the Student single tailed t test value for the 99% confidence level with

n-1 degrees of freedom and S denotes the standard deviation obtained from n replicate analyses. For the n=7 replicates used to determine the MDL,  $t=3.143$ . The results for the standard deviations and MDL are in Appendix 1.

## 12.2 Reporting Limit (RL)

Reporting limit (RL) refers to a level at which reliable quantitative results may be obtained. The MDL is used as a guide to determine the RL. The RL is chosen in a range 1-5 times the MDL, as per client requirement. The reporting limits for the selected rice herbicides are 0.05ppb.

## 12.3 Method Validation

The method validation consisted of 5 sample sets in background matrix. Each set included five levels of fortification for the selected rice herbicides (0.10, 0.25, 0.5, 1.25 and 2.5 ppb) and a method blank. All spikes and method blanks are processed through the entire analytical method. Recoveries for the selected rice herbicides are tabulated in Appendix 2.

## 12.4 Control Charts and Limits

Control charts were generated using the data from the method validation. The upper and lower warning and control limits are set at  $\pm 2$  and 3 standard deviations of the average percent recovery, respectively, shown in Appendix 2.

## 12.5 Acceptance Criteria

12.5.1 Each set of samples will have a matrix blank and a spiked matrix sample.

12.5.2 For positive results the retention time shall not vary from standards more than  $\pm 0.1$  minute.

12.5.3 Presence of both Qual and Quan ion.

12.5.4 The recoveries of the matrix spikes shall be within the control limits. See Appendix 2.

12.5.5 The sample shall be diluted if result exceed 10% of the highest calibration standard on the curve.



12.5.6 The relative abundances of structurally significant ions used for confirmation must be within  $\pm 30\%$  relative when compared to a standard injected during the same run.

13. Calculations:

Quantitation is based on an external standard (ESTD) calculation using either the peak area or height. The software uses linear or quadratic curve fit. Alternatively, at the chemist's discretion, concentrations may be calculated using the response factor for the standard whose value is  $< 30\%$  to the level in the sample.

$$\text{ppb} = \frac{(\text{sample peak area or ht}) \times (\text{std conc}) \times (\text{std vol. injected}) \times (\text{final vol of sample}) (1000 \mu\text{L/mL})}{(\text{std. peak area or ht}) \times (\text{sample vol injected}) \times (\text{sample wt (g)})}$$

14. Reporting Procedure:

Sample results are reported out according to the client's analytical laboratory specification sheets.

15. Discussion:

Diclofenac-d4 was used as an internal standard added at the end of extraction to calculate results. Both external standard (ESTD) and internal standard (ISTD) were used to calculate results during this method validation. We found that the quantitation based on an external standard (ESTD) with matrix matching standards worked better.

16. References:

- 16.1 Tsai, Cindy, *Determination of Selected Rice Herbicides in Water Samples by Solid Phase Extraction (SPE) and LC/MSMS*, Department of Fish and Wildlife, SOP: WPCL-LC-009
- 16.2 Fitch, P., Tran, D., updated by Lee, P., Hsu, J., White, J. *Determination of Atrazine, Bromacil, Cyanazine, Diuron, Hexazinone, Metribuzin, Norflurazon, Prometon, Prometryn, Simazine, Deethyl Atrazine (DEA), Deisopropyl Atrazine (ACET), Diamino Chlorotrazine ( DACT), Tebuthiuron and the metabolites Tebuthiuron-104, Tebuthiuron-106, Tebuthiuron-107 and Tebuthiuron-108 in Well and River Water By Liquid Chromatography- triple quadrupole mass spectrometry (LC/MS/MS)*. 2014, Environmental Monitoring Method, Center for Analytical Chemistry, CDFA

## Appendix 1

### The Determination of Method Detection Limit (MDL) and Reporting Limit (RL) in Groundwater

calculated without internal std spike 0.1ppb

| Compound Name       | Spk 1 | Spk 2  | Spk 3  | Spk 4  | Spk 5  | Spk 6  | Spk 7  | SD       | MDL     |
|---------------------|-------|--------|--------|--------|--------|--------|--------|----------|---------|
| Bensulfuron-methyl  | 0.117 | 0.103  | 0.104  | 0.0957 | 0.0975 | 0.100  | 0.0997 | 0.007048 | 0.02215 |
| Bispyribac-sodium   | 0.115 | 0.0988 | 0.101  | 0.0925 | 0.0960 | 0.0972 | 0.0950 | 0.007412 | 0.02330 |
| Clomazone           | 0.115 | 0.101  | 0.102  | 0.0936 | 0.0981 | 0.0969 | 0.0983 | 0.006876 | 0.02161 |
| Halosulfuron-methyl | 0.115 | 0.0976 | 0.0999 | 0.0903 | 0.0948 | 0.0941 | 0.0959 | 0.007974 | 0.02506 |
| Molinate            | 0.114 | 0.0974 | 0.102  | 0.0873 | 0.0966 | 0.0956 | 0.0976 | 0.008081 | 0.02540 |
| Orthosulfamuron     | 0.112 | 0.102  | 0.102  | 0.0934 | 0.0985 | 0.101  | 0.0977 | 0.005743 | 0.01805 |
| Penoxsulam          | 0.117 | 0.102  | 0.103  | 0.0966 | 0.0994 | 0.101  | 0.100  | 0.006622 | 0.02081 |
| Propanil            | 0.118 | 0.102  | 0.102  | 0.0932 | 0.101  | 0.102  | 0.104  | 0.007403 | 0.02327 |
| Propiconazole       | 0.117 | 0.0984 | 0.102  | 0.0934 | 0.0961 | 0.0994 | 0.0993 | 0.007645 | 0.02403 |
| Thiobencarb         | 0.114 | 0.0926 | 0.0951 | 0.0799 | 0.0898 | 0.0870 | 0.0943 | 0.010525 | 0.03308 |
| Simazine-d10        | 0.121 | 0.104  | 0.104  | 0.0972 | 0.102  | 0.104  | 0.102  | 0.007503 | 0.02358 |
| MCPA                | 0.116 | 0.102  | 0.101  | 0.0911 | 0.100  | 0.0994 | 0.100  | 0.007397 | 0.02325 |
| Triclopyr           | 0.117 | 0.102  | 0.0996 | 0.0890 | 0.0975 | 0.0980 | 0.100  | 0.008397 | 0.02639 |
| 2,4,5 T             | 0.116 | 0.0952 | 0.0982 | 0.0850 | 0.0965 | 0.0940 | 0.0982 | 0.009291 | 0.02920 |

Reporting Limit =0.05ppb for all compounds

**Appendix 2**

| Compound            | Set# | Spike Level |         |         |         |      | Control Limits |      |
|---------------------|------|-------------|---------|---------|---------|------|----------------|------|
|                     |      | 0.1 ppb     | 0.25ppb | 0.5 ppb | 1.25ppb | 5ppb | %              |      |
| Bensulfuron-methyl  | 1    | 95.5        | 97.2    | 91.4    | 92.8    | 89.6 | Mean:          | 95.5 |
|                     | 2    | 89.5        | 97.2    | 86.4    | 98.4    | 78.8 | SD:            | 7.95 |
|                     | 3    | 106         | 109     | 103     | 104     | 99.2 | UCL:           | 119  |
|                     | 4    | 92.9        | 89.6    | 76.2    | 96.8    | 96.0 | uwl            | 111  |
|                     | 5    | 98.8        | 91.2    | 105     | 103     | 100  | lwl            | 79.6 |
|                     |      |             |         |         |         |      | LCL:           | 71.6 |
| Bispyribac sodium   | 1    | 93.6        | 91.6    | 90.6    | 89.6    | 88.0 | Mean:          | 92.0 |
|                     | 2    | 85.7        | 92.0    | 83.2    | 95.2    | 74.0 | SD:            | 6.96 |
|                     | 3    | 101         | 100     | 98.0    | 100     | 95.2 | UCL:           | 113  |
|                     | 4    | 89.6        | 87.2    | 78.2    | 95.2    | 92.8 | uwl            | 106  |
|                     | 5    | 95.1        | 86.4    | 99.4    | 99.2    | 99.2 | lwl            | 78.1 |
|                     |      |             |         |         |         |      | LCL:           | 68.8 |
| Clomazone           | 1    | 94.5        | 93.6    | 90.4    | 87.2    | 88.4 | Mean:          | 92.5 |
|                     | 2    | 88.2        | 94.8    | 83.2    | 95.2    | 77.2 | SD:            | 7.74 |
|                     | 3    | 103         | 99.6    | 97.6    | 101     | 98.4 | UCL:           | 116  |
|                     | 4    | 89.0        | 84.0    | 74.6    | 95.2    | 93.6 | uwl            | 108  |
|                     | 5    | 96.9        | 83.2    | 101     | 102     | 100  | lwl            | 77.0 |
|                     |      |             |         |         |         |      | LCL:           | 69.3 |
| Halosulfuron methyl | 1    | 90.0        | 89.2    | 88.6    | 89.6    | 86.8 | Mean:          | 92.0 |
|                     | 2    | 85.6        | 95.2    | 84.8    | 95.2    | 74.8 | SD:            | 7.56 |
|                     | 3    | 101         | 101     | 100     | 102     | 96.8 | UCL:           | 115  |
|                     | 4    | 89.3        | 85.2    | 77.6    | 95.2    | 92.4 | uwl            | 107  |
|                     | 5    | 93.3        | 84.4    | 100     | 102     | 100  | lwl            | 76.9 |
|                     |      |             |         |         |         |      | LCL:           | 69.3 |
| Molinate            | 1    | 92.1        | 93.2    | 90.0    | 85.6    | 89.2 | Mean:          | 92.0 |
|                     | 2    | 86.3        | 96.4    | 86.6    | 95.2    | 76.8 | SD:            | 7.26 |
|                     | 3    | 99.9        | 96.4    | 97.2    | 100     | 95.2 | UCL:           | 114  |
|                     | 4    | 91.2        | 88.0    | 78.6    | 99.2    | 95.2 | uwl            | 106  |
|                     |      | 94.3        | 75.2    | 101     | 98.4    | 98.0 | lwl            | 77.5 |
|                     | 5    |             |         |         |         |      | LCL:           | 70.2 |

|                 |   |      |      |      |      |      |  |
|-----------------|---|------|------|------|------|------|--|
| Orthosulfamuron | 1 | 83.2 | 87.2 | 87.0 | 88.0 | 86.0 | Mean: 93.3<br>SD: 7.64<br>UCL: 116<br>uwl: 109<br>lwl: 78.1<br>LCL: 70.4 |
|                 | 2 | 93.6 | 108  | 92.8 | 104  | 86.4 |  |
|                 | 3 | 95.0 | 103  | 101  | 100  | 98.4 |  |
|                 | 4 | 86.1 | 85.6 | 77.4 | 97.6 | 94.0 |  |
|                 | 5 | 89.4 | 90.4 | 99.0 | 99.2 | 101  |  |
|                 |   |      |      |      |      |      |  |
| Penoxsulam      | 1 | 92.2 | 90.0 | 88.6 | 88.6 | 87.6 | Mean: 93.2<br>SD: 6.99<br>UCL: 114<br>uwl: 107<br>lwl: 79.2<br>LCL: 72.2 |
|                 | 2 | 92.7 | 93.6 | 84.2 | 95.2 | 80.0 |  |
|                 | 3 | 104  | 99.6 | 96.8 | 100  | 97.2 |  |
|                 | 4 | 93.9 | 85.6 | 76.8 | 97.6 | 94.8 |  |
|                 | 5 | 98.4 | 88.8 | 101  | 102  | 101  |  |
|                 |   |      |      |      |      |      |  |
| Propanil        | 1 | 90.4 | 91.6 | 87.2 | 88.0 | 87.6 | Mean: 91.6<br>SD: 7.97<br>UCL: 116<br>uwl: 108<br>lwl: 75.7<br>LCL: 67.7 |
|                 | 2 | 81.9 | 94.8 | 84.6 | 91.2 | 72.0 |  |
|                 | 3 | 101  | 101  | 93.8 | 100  | 95.6 |  |
|                 | 4 | 88.6 | 86.0 | 76.6 | 95.2 | 92.4 |  |
|                 | 5 | 96.7 | 86.4 | 101  | 101  | 102  |  |
|                 |   |      |      |      |      |      |  |
| Propiconazole   | 1 | 92.3 | 90.4 | 89.0 | 90.4 | 88.8 | Mean: 92.5<br>SD: 7.91<br>UCL: 116<br>uwl: 108<br>lwl: 76.6<br>LCL: 68.7 |
|                 | 2 | 88.3 | 94.4 | 83.8 | 97.6 | 76.8 |  |
|                 | 3 | 104  | 102  | 98.0 | 102  | 96.0 |  |
|                 | 4 | 93.6 | 87.6 | 78.0 | 96.0 | 80.4 |  |
|                 | 5 | 95.9 | 82.4 | 101  | 101  | 102  |  |
|                 |   |      |      |      |      |      |  |
| Thiobencarb     | 1 | 89.1 | 88.0 | 85.8 | 84.8 | 88.0 | Mean: 88.2<br>SD: 8.95<br>UCL: 115<br>uwl: 106<br>lwl: 70.3<br>LCL: 61.3 |
|                 | 2 | 73.7 | 92.4 | 81.2 | 81.6 | 70.0 |  |
|                 | 3 | 105  | 98.8 | 98.4 | 97.6 | 93.6 |  |
|                 | 4 | 83.8 | 82.8 | 73.8 | 91.2 | 84.8 |  |
|                 | 5 | 89.1 | 78.4 | 94.4 | 98.4 | 99.2 |  |
|                 |   |      |      |      |      |      |  |
| Simazine-d10    | 1 | 91.5 | 83.6 | 80.8 | 76.7 | 78.4 | Mean: 89.2<br>SD: 8.64<br>UCL: 115<br>uwl: 106<br>lwl: 71.9<br>LCL: 63.2 |
|                 | 2 | 87.1 | 94.0 | 88.0 | 99.2 | 77.2 |  |
|                 | 3 | 104  | 94.8 | 88.4 | 93.6 | 93.6 |  |
|                 | 4 | 93.9 | 80.8 | 74.8 | 92.8 | 89.2 |  |
|                 | 5 | 98.3 | 76.4 | 92.8 | 94.4 | 105  |  |
|                 |   |      |      |      |      |      |  |

|           |   |      |      |      |      |      |       |      |
|-----------|---|------|------|------|------|------|-------|------|
| MCPA      | 1 | 96.7 | 93.6 | 91.0 | 92.8 | 91.6 | Mean: | 93.5 |
|           | 2 | 85.2 | 95.2 | 84.2 | 92.0 | 72.0 | SD:   | 7.63 |
|           | 3 | 105  | 102  | 97.6 | 100  | 98.4 | UCL:  | 116  |
|           | 4 | 93.8 | 88.0 | 79.2 | 96.8 | 93.6 | uwl   | 109  |
|           | 5 | 98.5 | 88.4 | 101  | 101  | 101  | lwl   | 78.3 |
|           |   |      |      |      |      |      | LCL:  | 70.7 |
| Triclopyr | 1 | 95.5 | 95.6 | 91.8 | 95.2 | 94.0 | Mean: | 94.8 |
|           | 2 | 83.4 | 93.2 | 83.6 | 92.0 | 71.2 | SD:   | 8.24 |
|           | 3 | 104  | 107  | 100  | 102  | 101  | UCL:  | 120  |
|           | 4 | 94.8 | 87.2 | 84.0 | 99.2 | 96.8 | uwl   | 111  |
|           | 5 | 97.6 | 90.4 | 104  | 102  | 104  | lwl   | 78.3 |
|           |   |      |      |      |      |      | LCL:  | 70.1 |
| 2,4,5,T   | 1 | 96.7 | 94.4 | 91.4 | 94.4 | 92.4 | Mean: | 93.3 |
|           | 2 | 77.6 | 92.8 | 84.0 | 86.4 | 69.2 | SD:   | 8.81 |
|           | 3 | 104  | 104  | 99.6 | 102  | 98.8 | UCL:  | 120  |
|           | 4 | 91.2 | 88.8 | 81.0 | 97.6 | 91.6 | uwl   | 111  |
|           | 5 | 97.9 | 88.4 | 102  | 102  | 104  | lwl   | 75.7 |
|           |   |      |      |      |      |      | LCL:  | 66.9 |

**Written By:**

Original signed by

1/31/2020

---

Jean Hsu  
Sr. Environmental Scientist (Specialist)

---

Date

**Written By:**

Original signed by

1/31/2020

---

Jane White  
Sr. Environmental Scientist (Specialist)

---

Date

**Approved By:**

Original signed by

1/31/2020

---

Stephen Siegel  
Sr. Environmental Scientist Supervisor

---

Date

**Approved By:**

Original signed by

1/31/2020

---

Maryam Khosravifard  
Environmental Program Manager I

---

Date

**Approved By:**

Original signed by

3/5/2020

---

Sarva Balachandra  
Quality Assurance Officer

---

Date

**Revision Log:**

