

Status: **Final** Approved by:

---

**Mui Koltunov, Ph.D., Branch Chief**  
**Analytical Chemistry Branch, Pasadena**

---

## **Pesticide Analysis by Liquid Chromatography Triple Quadrupole Mass Spectrometry (LC-MS-MS)**

---

Written by: **Jessica Dawlaty, PhD**

### **SCOPE AND APPLICABILITY**

Monitoring pesticides in wastewater, surface water, and drinking water is an important process of water quality control. Pesticides may be harmful to humans, animals, and the environment because of their toxicity. These pesticides may be used in crops, buildings, ornamental plants, and lawns. The routes of exposure include oral, inhalation, and dermal. To minimize the risk of exposure to these chemicals, Congress passed the Food Quality Protection Act of 1996 to allow EPA to regulate the use of pesticides.

This technical system procedure (TSP) describes steps to determine the concentration of a group of mainly pyrethroids pesticides from extracts prepared from aqueous liquid matrix. Target pesticides in wastewater samples are extracted using TSP-04.3535.00, Solid-Phase Extraction for Pesticides, with target pesticides listed in Table 1. The extracts are then divided into two aliquots and analyzed by gas chromatography quadrupole time-of-flight mass spectrometry (GC-QToF-MS) and by liquid chromatography triple quadrupole mass spectrometry (LC-MS-MS). This procedure focuses on data acquisition, data analysis and reporting for LC-MS-MS.

### **PRINCIPLE**

This method quantifies pesticides using liquid chromatography triple quadrupole mass spectrometry (LC-MS-MS). The compounds of interest are separated using high performance liquid chromatography (HPLC) on a reverse phase C-18 column. After separation, the pesticides are ionized by electrospray ionization and directed into the mass spectrometer. The triple quadrupole mass spectrometer is operated in dynamic multiple reaction monitoring (dMRM) mode. In dMRM, each compound has a retention time specific time window for acquisition. During acquisition, a precursor ion mass is selected in the first quadrupole, fragmented in the collision cell, compound specific product ions are selected in the last quadrupole, and finally sent to the detector.

### **SAFETY**

#### **Precautions**

Ethyl acetate, methanol, and ammonium fluoride are hazardous and volatile, avoid direct inhalation of the solvents. Be aware of the location of exit door(s), eye wash and shower station(s), spill kit, first aid kit and closest phone in case of emergency.

### **Personal Protective Equipment (PPE)**

Wear lab coat, glasses, and gloves at all times. Fume hood is used whenever possible.

### **INTERFERENCES**

Ensure glassware and consumables are free of contamination. Plastics may cause phthalate contamination.

### **PRESERVATION AND HOLDING TIMES**

#### **Storage and Holding Times**

Aqueous liquid samples are stored at <6°C in Teflon-sealed glass containers. Aqueous liquid samples have 14 days to be extracted and extracts should be analyzed within 40 days.

### **EQUIPMENT AND SUPPLIES**

#### **Hardware**

Agilent HPLC LC-MS-MS System: Agilent 1290 Infinity II HPLC, Agilent 6470 Triple Quadrupole Mass Spectrometer, Agilent MassHunter Software.  
Zorbax RRHD Eclipse Plus C18 column: 2.1 × 100 mm, 1.8 μm

#### **Instrumentation Conditions**

LC-MS-MS dMRM Transitions for Pesticide Targets, Surrogates, and Internal Standards, Table 1  
Acquisition Parameters for LC-QQQ, Table 2

#### **Consumables**

Autosampler vials – 2 mL with caps with PTFE-lined septa, Nitrile gloves. Avoid vinyl gloves. Gas tight syringes: 10 μL, 25 μL, 100 μL, 250 μL and 1 mL.

### **REFERENCE MATERIALS AND REAGENTS**

#### **Reference Materials (RM)**

Custom pesticide standard from Accustandard. See Table 3.  
Refer to FRM-07.0329.00, Pesticide Reference Materials (RMs) Preparation Log, for list of RMs.  
Refer to FRM-07.0330.00.00, Pesticide Working Reference Materials (WRMs) Preparation Log, for preparation and list of WRMs.  
Refer to FRM-07.0328.00, Pesticide Calibration and Verification Solution Preparation for LC-QQQ, for preparation and list of calibration and verification WRMs

#### **Reagents**

Methanol, HPLC grade, or equivalent  
Water, HPLC grade, or equivalent

Isopropanol, HPLC grade, or equivalent  
Ammonium fluoride, mass spectrometry grade, or equivalent

## **METHOD PROCEDURE**

### **Instrument Requirements**

Refer to QSP-02.0048.00, Instrument Checks, for requirements.

### **Batch Quality Control (QC) Requirements**

Refer to QSP-02.0024.00, Quality Control, for requirements

### **Sample Preparation**

Samples are extracted using TSP-04.3535.00, Automated Solid-Phase Extraction for Pesticides, Allow sample extracts to warm to room temperature, if necessary. Prior to analysis, add 10 µL of 5 µg/mL Internal Standard Working Solution to each 0.5 mL extract. Place the samples in the autosampler.

### **MassHunter Data Acquisition and Analysis**

Refer to TSP-06.0007.00, MassHunter Data Acquisition

### **Data Reporting**

Complete the applicable report: templates are located in the “U” drive/ECL Pasadena/Report Templates and Archive/Report Templates. Refer to TSP-06.0007.00, MassHunter Data Acquisition Section 4.5. for a list of reports from MassHunter that are required. Make copies of all supporting records. Refer to Report Package Contents (check-off) on FRM-07.0106.00, Initial Level Review Check List, when compiling the report package.

## **MAINTENANCE AND TROUBLESHOOTING**

### **Mass Spectrometer Maintenance**

Prior to sample analysis, rinse mass spectrometer source with 50% isopropanol in water, followed by methanol. After analysis, clean the nebulizer. Clean the mass spectrometer source as needed per manufacturer’s instructions.

### **Waste Disposal**

Solvent waste is poured into the appropriate waste bottle. Extracts should be stored in the sample refrigerator for a minimum of 40 days before being disposed of in the solvent waste bottle. Glass reference material vials are dried in the fume hood and disposed in a waste bag.

## **REFERENCES**

United States Environmental Protection Agency, “*Method 3500C - Organic Extraction and Sample Preparation*”, SW-846, 2007.

California Environmental Protection Agency  
Department of Toxic Substances Control  
Environmental Chemistry Laboratory  
UNCONTROLLED DOCUMENT - DO NOT COPY

DCN: 06.0016.00  
Revision No: 0  
August 5, 2022  
Page 4 of 10

United States Environmental Protection Agency, "*Method 8000C - Determinative Chromatographic Separations*", SW-846, 2003.

**Table 1. LC-MS-MS dMRM Transitions for Pesticide Targets, Surrogates, and Internal Standards**

| Target Analyte          | Retention Time (min) | Retention Window (min) | Precursor Ion (m/z) | Product Ion (m/z) | Fragmentor (V) | Collision Energy (V) | Cell Acceleration (V) | Polarity |
|-------------------------|----------------------|------------------------|---------------------|-------------------|----------------|----------------------|-----------------------|----------|
| Bifenthrin              | 9                    | 0.6                    | 442                 | 181               | 100            | 16                   | 2                     | Positive |
| Bifenthrin              | 9                    | 0.6                    | 442                 | 166               | 100            | 35                   | 2                     | Positive |
| Bifenthrin-D6 (ISTD)    | 9                    | 0.6                    | 446.2               | 181.1             | 68             | 16                   | 2                     | Positive |
| Bifenthrin-D6 (ISTD)    | 9                    | 0.6                    | 446.2               | 166.1             | 68             | 50                   | 2                     | Positive |
| Bioallethrin            | 5.26                 | 0.4                    | 303.2               | 135               | 80             | 10                   | 3                     | Positive |
| Bioallethrin            | 5.26                 | 0.4                    | 303.2               | 107               | 80             | 20                   | 3                     | Positive |
| Chlorpyrifos            | 5.45                 | 0.4                    | 349.9               | 197.92            | 100            | 20                   | 5                     | Positive |
| Chlorpyrifos            | 5.45                 | 0.4                    | 349.9               | 96.95             | 100            | 41                   | 5                     | Positive |
| Chlorpyrifos-D10 (Surr) | 5.45                 | 0.4                    | 360                 | 198.9             | 102            | 20                   | 2                     | Positive |
| Chlorpyrifos-D10 (Surr) | 5.45                 | 0.4                    | 360                 | 99                | 102            | 44                   | 2                     | Positive |
| Cyfluthrin              | 6.4                  | 0.8                    | 451.3               | 434.1             | 58             | 4                    | 2                     | Positive |
| Cyfluthrin              | 6.4                  | 0.8                    | 451.3               | 191               | 58             | 13                   | 2                     | Positive |
| Cyfluthrin-13C6 (ISTD)  | 6.68                 | 0.8                    | 439.1               | 422.1             | 68             | 4                    | 2                     | Positive |
| Cyfluthrin-13C6 (ISTD)  | 6.68                 | 0.8                    | 439.1               | 197               | 68             | 16                   | 2                     | Positive |
| Cyhalothrin             | 6.68                 | 0.8                    | 467.13              | 450               | 100            | 5                    | 3                     | Positive |
| Cyhalothrin             | 6.68                 | 0.8                    | 467.13              | 225               | 100            | 10                   | 3                     | Positive |
| Cyhalothrin-D6 (ISTD)   | 6.51                 | 0.8                    | 473.2               | 456               | 100            | 5                    | 2                     | Positive |
| Cyhalothrin-D6 (ISTD)   | 6.51                 | 0.8                    | 473.2               | 231               | 58             | 16                   | 2                     | Positive |
| Cypermethrin            | 6.67                 | 0.8                    | 433.3               | 416               | 90             | 7                    | 2                     | Positive |
| Cypermethrin            | 6.67                 | 0.8                    | 433.3               | 191               | 90             | 16                   | 2                     | Positive |
| Cypermethrin-D6 (ISTD)  | 6.62                 | 0.8                    | 439.1               | 422.1             | 53             | 8                    | 2                     | Positive |
| Cypermethrin-D6 (ISTD)  | 6.62                 | 0.8                    | 439.1               | 197               | 53             | 16                   | 2                     | Positive |
| Deltamethrin            | 7                    | 0.8                    | 523                 | 280.9             | 70             | 15                   | 3                     | Positive |
| Deltamethrin            | 7                    | 0.8                    | 523                 | 181               | 70             | 50                   | 3                     | Positive |
| Deltamethrin-D6 (ISTD)  | 6.8                  | 0.6                    | 527                 | 510               | 53             | 8                    | 2                     | Positive |
| Deltamethrin-D6 (ISTD)  | 6.8                  | 0.6                    | 527                 | 284.9             | 53             | 16                   | 2                     | Positive |
| Diuron-D6 (Surr)        | 3.02                 | 0.4                    | 239.1               | 78.1              | 81             | 32                   | 3                     | Positive |
| Diuron-D6 (Surr)        | 3.02                 | 0.4                    | 239.1               | 52.1              | 81             | 20                   | 6                     | Positive |
| Esfenvalerate           | 7.12                 | 0.5                    | 437.16              | 167               | 100            | 10                   | 3                     | Positive |
| Esfenvalerate           | 7.12                 | 0.5                    | 437.16              | 125               | 100            | 25                   | 3                     | Positive |
| Etofenprox              | 8.5                  | 0.6                    | 394.2               | 177.3             | 90             | 8                    | 2                     | Positive |
| Etofenprox              | 8.5                  | 0.6                    | 394.2               | 107.1             | 90             | 40                   | 2                     | Positive |
| Etofenprox-D5 (ISTD)    | 8.45                 | 0.6                    | 399.3               | 182.1             | 53             | 12                   | 2                     | Positive |
| Etofenprox-D5 (ISTD)    | 8.45                 | 0.6                    | 399.3               | 108               | 53             | 50                   | 2                     | Positive |
| Fenpropathrin           | 6.2                  | 0.4                    | 350.2               | 125.1             | 115            | 10                   | 3                     | Positive |
| Fenpropathrin           | 6.2                  | 0.4                    | 350.2               | 97.1              | 115            | 32                   | 3                     | Positive |
| Fenpropathrin-D6 (ISTD) | 6.2                  | 0.4                    | 373.2               | 356.2             | 53             | 4                    | 2                     | Positive |
| Fenpropathrin-D6 (ISTD) | 6.2                  | 0.4                    | 373.2               | 131.1             | 53             | 16                   | 2                     | Positive |
| Fipronil                | 3.6                  | 0.4                    | 434.9               | 330               | 106            | 16                   | 2                     | Negative |

| Target Analyte                          | Retention Time (min) | Retention Window (min) | Precursor Ion (m/z) | Product Ion (m/z) | Fragmentor (V) | Collision Energy (V) | Cell Acceleration (V) | Polarity |
|---|----------------------|------------------------|---------------------|-------------------|----------------|----------------------|-----------------------|----------|
| Fipronil                                | 3.6                  | 0.4                    | 434.9               | 250               | 106            | 32                   | 4                     | Negative |
| Fipronil-[13C]4[15N]2 (ISTD)            | 3.6                  | 0.4                    | 440.9               | 336               | 115            | 16                   | 2                     | Negative |
| Fipronil-[13C]4[15N]2 (ISTD)            | 3.6                  | 0.4                    | 440.9               | 252               | 115            | 32                   | 4                     | Negative |
| Fipronil-amide                          | 3                    | 0.5                    | 454.96              | 385.9             | 121            | 20                   | 2                     | Negative |
| Fipronil-amide                          | 3                    | 0.5                    | 454.96              | 290               | 121            | 40                   | 2                     | Negative |
| Fipronil-desulfinyl                     | 3.48                 | 0.4                    | 387                 | 351               | 86             | 12                   | 2                     | Negative |
| Fipronil-desulfinyl                     | 3.48                 | 0.4                    | 387                 | 282               | 86             | 36                   | 2                     | Negative |
| Fipronil-desulfinyl amide               | 2.85                 | 0.4                    | 405                 | 369               | 82             | 8                    | 2                     | Negative |
| Fipronil-desulfinyl amide               | 2.85                 | 0.4                    | 405                 | 329               | 82             | 20                   | 2                     | Negative |
| Fipronil-desulfinyl-[13C]4[15N]2 (ISTD) | 3.48                 | 0.4                    | 393                 | 357               | 101            | 12                   | 2                     | Negative |
| Fipronil-desulfinyl-[13C]4[15N]2 (ISTD) | 3.48                 | 0.4                    | 393                 | 288               | 101            | 36                   | 2                     | Negative |
| Fipronil-sulfide                        | 3.66                 | 0.4                    | 418.9               | 383               | 101            | 12                   | 2                     | Negative |
| Fipronil-sulfide                        | 3.66                 | 0.4                    | 418.9               | 262               | 101            | 32                   | 2                     | Negative |
| Fipronil-sulfide-[13C]4[15N]2 (ISTD)    | 3.66                 | 0.4                    | 424.9               | 389               | 115            | 12                   | 2                     | Negative |
| Fipronil-sulfide-[13C]4[15N]2 (ISTD)    | 3.66                 | 0.4                    | 424.9               | 265               | 115            | 32                   | 2                     | Negative |
| Fipronil-sulfone                        | 3.78                 | 0.4                    | 450.9               | 415               | 115            | 16                   | 2                     | Negative |
| Fipronil-sulfone                        | 3.78                 | 0.4                    | 450.9               | 282               | 115            | 28                   | 4                     | Negative |
| Fipronil-sulfone-[13C]4[15N]2 (ISTD)    | 3.78                 | 0.4                    | 456.9               | 421               | 130            | 16                   | 2                     | Negative |
| Fipronil-sulfone-[13C]4[15N]2 (ISTD)    | 3.78                 | 0.4                    | 456.9               | 288               | 130            | 32                   | 4                     | Negative |
| Fonofos (Surr)                          | 3.99                 | 0.4                    | 247.04              | 137.02            | 80             | 5                    | 3                     | Positive |
| Fonofos (Surr)                          | 3.99                 | 0.4                    | 247.04              | 108.99            | 80             | 15                   | 3                     | Positive |
| Imidacloprid                            | 1.78                 | 0.8                    | 256                 | 208.9             | 80             | 12                   | 4                     | Positive |
| Imidacloprid                            | 1.78                 | 0.8                    | 256                 | 175               | 80             | 12                   | 4                     | Positive |
| Imidacloprid-D4 (ISTD)                  | 1.75                 | 0.8                    | 260.1               | 213               | 67             | 16                   | 2                     | Positive |
| Imidacloprid-D4 (ISTD)                  | 1.75                 | 0.8                    | 260.1               | 179.1             | 67             | 24                   | 2                     | Positive |
| Linuron-D6 (ISTD)                       | 3.2                  | 0.4                    | 255.1               | 185               | 91             | 20                   | 3                     | Positive |
| Linuron-D6 (ISTD)                       | 3.2                  | 0.4                    | 255.1               | 159.9             | 91             | 20                   | 3                     | Positive |
| Methomyl-D3 (Surr)                      | 1.2                  | 0.6                    | 166.1               | 105.9             | 48             | 8                    | 3                     | Positive |
| Methomyl-D3 (Surr)                      | 1.2                  | 0.6                    | 166.1               | 88                | 48             | 8                    | 3                     | Positive |
| Novaluron                               | 4.52                 | 0.4                    | 493                 | 158.1             | 90             | 16                   | 3                     | Positive |
| Novaluron                               | 4.52                 | 0.4                    | 493                 | 141.1             | 90             | 56                   | 3                     | Positive |
| Permethrin                              | 7.71                 | 0.4                    | 391.09              | 355               | 100            | 5                    | 3                     | Positive |
| Permethrin                              | 7.71                 | 0.4                    | 391.09              | 183               | 100            | 5                    | 3                     | Positive |
| Permethrin-D6 (ISTD)                    | 7.66                 | 0.4                    | 414.1               | 361.1             | 48             | 4                    | 2                     | Positive |

| Target Analyte                | Retention Time (min) | Retention Window (min) | Precursor Ion (m/z) | Product Ion (m/z) | Fragmentor (V) | Collision Energy (V) | Cell Acceleration (V) | Polarity |
|-------------------------------|----------------------|------------------------|---------------------|-------------------|----------------|----------------------|-----------------------|----------|
| Permethrin-D6 (ISTD)          | 7.66                 | 0.4                    | 414.1               | 183               | 48             | 20                   | 2                     | Positive |
| Prallethrin                   | 4.4                  | 0.4                    | 301.18              | 169               | 80             | 5                    | 3                     | Positive |
| Prallethrin                   | 4.4                  | 0.4                    | 301.18              | 105               | 80             | 20                   | 3                     | Positive |
| Propoxur                      | 2.8                  | 0.4                    | 210.11              | 168.1             | 48             | 4                    | 4                     | Positive |
| Propoxur                      | 2.8                  | 0.4                    | 210.11              | 153.1             | 48             | 4                    | 3                     | Positive |
| Propoxur-D3 (ISTD)            | 2.8                  | 0.4                    | 213.1               | 171.1             | 106            | 4                    | 4                     | Positive |
| Propoxur-D3 (ISTD)            | 2.8                  | 0.4                    | 213.1               | 111               | 106            | 16                   | 3                     | Positive |
| Pyrethrin 1                   | 6.21                 | 0.4                    | 329.2               | 161               | 67             | 8                    | 3                     | Positive |
| Pyrethrin 1                   | 6.21                 | 0.4                    | 329.2               | 133.1             | 67             | 20                   | 3                     | Positive |
| Pyriproxyfen                  | 5.3                  | 0.4                    | 322.2               | 185               | 110            | 20                   | 2                     | Positive |
| Pyriproxyfen                  | 5.3                  | 0.4                    | 322.2               | 96                | 110            | 12                   | 2                     | Positive |
| Simazine-D5 (Surr)            | 2.8                  | 0.4                    | 207.1               | 137               | 105            | 20                   | 3                     | Positive |
| Simazine-D5 (Surr)            | 2.8                  | 0.4                    | 207.1               | 129               | 105            | 20                   | 4                     | Positive |
| Tau-Fluvalinate               | 7.8                  | 0.4                    | 503.13              | 208.08            | 80             | 15                   | 3                     | Positive |
| Tau-Fluvalinate               | 7.8                  | 0.4                    | 503.13              | 181.06            | 80             | 25                   | 3                     | Positive |
| Tetrachlorvinphos             | 3.75                 | 0.4                    | 364.9               | 203.9             | 120            | 40                   | 3                     | Positive |
| Tetrachlorvinphos             | 3.75                 | 0.4                    | 364.9               | 127               | 120            | 16                   | 3                     | Positive |
| Tetrachlorvinphos-13C6 (ISTD) | 3.75                 | 0.4                    | 370.9               | 209.9             | 121            | 50                   | 2                     | Positive |
| Tetrachlorvinphos-13C6 (ISTD) | 3.75                 | 0.4                    | 370.9               | 127               | 121            | 16                   | 2                     | Positive |
| Tetramethrin                  | 4.9                  | 0.8                    | 332.19              | 164.07            | 100            | 15                   | 3                     | Positive |
| Tetramethrin                  | 4.9                  | 0.8                    | 332.19              | 135               | 100            | 15                   | 3                     | Positive |

**Table 2. Acquisition Parameters for LC-MS-MS**

| <b>LC Parameters</b>            |                               |
|---------------------------------|-------------------------------|
| Injection volume                | 2 $\mu$ L                     |
| Needle wash                     | Flush port, 6 s               |
| Column temperature              | 55 $^{\circ}$ C               |
| Maximum pressure                | 800 bar                       |
| Flow rate                       | 0.4 mL/min                    |
| Eluent A                        | Water, 1 mM ammonium fluoride |
| Eluent B                        | Methanol                      |
| <b>Gradient</b>                 |                               |
| Time                            | % B                           |
| 0.00                            | 30 %                          |
| 1.00                            | 30 %                          |
| 2.00                            | 75 %                          |
| 10.00                           | 90 %                          |
| 10.20                           | 100 %                         |
| Stop time                       | 12.00 min                     |
| Post Time                       | 3.00 min                      |
| <b>Source Parameters</b>        |                               |
| Ion Source                      | AJS ESI                       |
| Gas Temp ( $^{\circ}$ C)        | 200                           |
| Gas Flow (L/min)                | 13                            |
| Nebulizer (psi)                 | 55                            |
| Sheath Gas Temp ( $^{\circ}$ C) | 200                           |
| Sheath Gas Flow (L/min)         | 10                            |
| Capillary (V)                   | 5000 (+) / 3500 (-)           |
| Nozzle (V)                      | 0                             |



**Table 3. Custom Pesticide Reference Materials**

| <b>Compound</b>         | <b>CAS ID</b> | <b>Concentration<br/>in mix (ug/mL)</b> |
|-------------------------|---------------|---|
| Bifenthrin              | 82657-04-3    | 100                                     |
| S-Bioallethrin          | 28434-00-6    | 50                                      |
| Chlorothalonil          | 1897-45-6     | 50                                      |
| Chlorpyrifos            | 2921-88-2     | 50                                      |
| Cyfluthrin              | 68359-37-5    | 50                                      |
| Cyhalothrin             | 68085-85-8    | 100                                     |
| Cypermethrin            | 52315-07-8    | 50                                      |
| Cyphenothrin            | 39515-40-7    | 250                                     |
| Deltamethrin            | 52918-63-5    | 250                                     |
| Esfenvalerate (Asana)   | 66230-04-4    | 100                                     |
| Etofenprox              | 80844-07-1    | 50                                      |
| Fenpropathrin (Danitol) | 39515-41-8    | 250                                     |
| Fipronil                | 120068-37-3   | 50                                      |
| Fipronil-desulfinyl     | 205650-65-3   | 50                                      |
| Fipronil-sulfide        | 120067-83-6   | 50                                      |
| Fipronil-sulfone        | 120068-36-2   | 50                                      |
| Imidacloprid            | 138261-41-3   | 50                                      |
| Novaluron               | 116714-46-6   | 50                                      |
| Permethrin              | 52645-53-1    | 2500                                    |
| Phenothrin              | 26002-80-2    | 2500                                    |
| Prallethrin             | 23031-36-9    | 50                                      |
| Propoxur (Baygon)       | 114-26-1      | 50                                      |
| Pyriproxyfen            | 95737-68-1    | 50                                      |
| Tau-fluvalinate         | 102851-06-9   | 50                                      |
| Tetrachlorvinphos       | 22248-79-9    | 50                                      |
| Tetramethrin            | 7696-12-0     | 50                                      |

**RE-APPROVAL**

**Signatures**

**Date**

---

---

---

---

---

---

---

---

---

---